

SPIROVIT SP-1 G2


Portable spirometer with ultrasonic technology
Innovation for optimal lung care



SPIROVIT SP-1 G2

SCHILLER's SPIROVIT SP-1 G2 is the all-in-one solution in lung function testing – portable/ desktop spirometry for large and small healthcare providers.

- Ultrasound technology: calibration-free, no warm up time, no moving parts.
- Real-time information on colour LCD



ULTRASOUND TECHNOLOGY

Thanks to GANSHORN's ultrasound technology, recordings are highly accurate and the sensor is maintenance- and calibration-free. Furthermore, the system is independent of gas composition, air pressure, temperature and humidity. The disposable ScoutTubes ensure maximum hygiene and patient safety.


DISPOSABLE AND SUSTAINABLE

To ensure a sustainable product, SCHILLER compensates the CO₂ emitted during production, packaging, shipment and disposal of the ScoutTubes.



REAL-TIME DISPLAY

The large colour LCD helps to ensure that the patient achieves a good result and the measurement is valid: information is displayed in real time, the spirometer provides feedback on the procedure's quality after each trial and an incentive screen is available for trials with children.



Incentive screen for paediatric spirometry

WIDE RANGE OF MEASUREMENTS INCLUDING INTERPRETATION

The spirometer offers FVC, SVC and MVV. Moreover, Pre and Post tests are possible. The SPIROVIT SP-1 G2 automatically interprets the results and provides information on a possible diagnosis. Results can be printed on site, using the integrated thermal printer.

AVAILABLE WITH TROLLEY


The SPIROVIT SP-1 G2 offers the major advantage of being able to perform spirometry wherever you are, without a PC. For added convenience, the system is also available with a trolley, which includes a bin for consumables such as ScoutTubes.

FVC, SVC, MVV, Pre/Post testing, wide range of predicted normal values

Export spirometry recordings from internal memory as PDF and store on external media

Colour LCD: information displayed in real time

High-resolution thermal printer



Ultrasonic sensor: calibration-free and maintenance-free

Highly hygienic sealed keyboard keeping dust, dirt and liquids out

Compliance with 2019 ATS/ERS Guidelines

Disposable ScoutTubes for maximum hygiene and patient safety





Americas
SCHILLER Americas Inc.
 Doral, Florida 33172

North America:
 Phone +1 786 845 06 20
 Fax +1 786 845 06 02
 sales@schilleramericas.com
 www.schilleramericas.com

Latin America & Caribbean
 Phone +1 305 591 11 21
 Fax +1 786 845 06 02
 sales@schilleramericas.com
 www.schilleramericas.com



Asia
SCHILLER Asia-Pacific / Malaysia
 52200 Kuala Lumpur, Malaysia
 Phone +603 6272 3033
 sales@schiller.com.my
 www.schiller-asia.com



Austria
SCHILLER Handelsgesellschaft m.b.H.
 A-4040 Linz
 Phone +43 732 70 99 0
 Fax +43 732 757 000
 sales@schiller.at
 www.schiller.at



China
Alfred Schiller (Beijing) Medical Equipment Co.,Ltd.
 100015 Beijing, China
 Phone +86 010 52007020
 Fax +86 010 52007020-8016
 info@schillerchina.com
 www.schillermedical.cn



Croatia
SCHILLER d.o.o.
 10000 Zagreb
 Phone +385 1 309 66 59
 Fax +385 1 309 66 60
 info@schillerzg.hr
 www.schiller.hr



France
SCHILLER Médical
 F-67160 Wissembourg
 Phone +33 3 88 63 36 00
 Fax +33 3 88 63 36 49
 info@schiller.fr
 www.schiller.fr



France (distribution France)
SCHILLER France S.A.S.
 F-77608 Bussy St Georges
 Phone +33 1 64 66 50 00
 Fax +33 1 64 66 50 10
 contact@schillerfrance.fr
 www.schiller-france.com



Germany
SCHILLER Medizintechnik GmbH
 D-85622 Feldkirchen b. München
 Phone +49 89 62 99 81 0
 Fax +49 89 609 50 90
 info@schillermed.de
 www.schillermed.de



India
SCHILLER Healthcare India Pvt. Ltd.
 Mumbai – 400 059, India
 Phone +91 22 6152 3333 / 2920 9141
 Fax +91 22 2920 9142
 sales@schillerindia.com
 www.schillerindia.com



Poland
SCHILLER Poland Sp. z o.o.
 PL-02-729 Warszawa
 Phone +48 22 843 20 89 / +48 22 647 35 90
 Fax +48 22 843 20 89
 schiller@schiller.pl
 www.schiller.pl



Russia & C.I.S.
AO SCHILLER.RU
 119049 Moscow, Russia
 Phone +7 (495) 970 11 33
 mail@schiller.ru
 www.schiller.ru



Serbia
SCHILLER d.o.o.
 11010 Beograd
 Phone +381 11 39 79 508
 Fax +381 11 39 79 518
 info@schiller.rs
 www.schiller.rs



Slovenia
SCHILLER d.o.o.
 2310 Slovenska Bistrica
 Phone +386 2 843 00 56
 Fax +386 2 843 00 57
 info@schiller.si
 www.schiller.si



Spain
SCHILLER ESPAÑA, S.A.
 ES-28232 Las Rozas/Madrid
 Phone +34 91 713 01 76
 Fax +34 91 355 79 33
 schiller@schiller.es
 www.schiller.es



Switzerland
SCHILLER-Reomed AG
 CH-8912 Obfelden
 Phone +41 44 744 30 00
 Fax +41 44 740 37 10
 sales@schiller-reomed.ch
 www.schiller-reomed.ch



Turkey
SCHILLER TÜRKIYE
 Okmeydani-Sisli – Istanbul
 Phone +90 212 210 8681 (pbx)
 Fax +90 212 210 8684
 info@schiller.com.tr
 www.schiller-turkiye.com



United Kingdom
SCHILLER UK Ltd.
 Anstruther, Fife KY10 3H
 Phone +44 16 190 72 565
 sales@schilleruk.com
 www.schilleruk.com

Device availability in your market is subject to regulatory approval.

